



White paper

Polarity-Extended Chromatography:

A holistic solution for the analysis of organic molecules in the aqueous environment


Stefan Bieber and Thomas Letzel

The gaining interest in persistent, mobile and toxic (PMT) and very persistent and very mobile (vPvM) substances in waterbodies resulted in the requirement to adjust the applied analytical techniques towards these compounds, which are often (very) polar, i.e. mostly with a $\log D(\text{pH}) \leq 0$ and some of them ionic. On the molecular level these (organic) compounds are perfectly suitable to analytical techniques like mass spectrometric detection and a prior chromatographic/electrophoretic separation.

This comment is a contribution to the special issue 'Persistent and Mobile Organic Compounds – An Environmental Challenge' of the Journal 'Analytical and Bioanalytical Chemistry' which we recommend to read [1]. Herewith, we would like to address and point towards chromatographic techniques which are available and useable for analyzing (simultaneously and combined) non-polar, polar and very polar molecules in one single run with so called 'polarity-extended chromatography'.



Commonly applied chromatographic techniques like gas chromatography and reversed phase liquid chromatography (RPLC) can be used for a chromatographic separation of non-polar and polar molecules. However, they reach their limits in analyzing very polar molecules. For such molecules, techniques like polar-modified RPLC, Hydrophilic Interaction Liquid Chromatography (HILIC), Ion Chromatography (IC), Capillary Electrophoresis (CE) and Supercritical Fluid Chromatography (SFC) can be applied to close the 'analytical gap' in the very polar region (see Figure below).



Polarity scheme for chromatographic (electrophoretic) separation techniques regarding the logD (logP) values of separable molecules and molecule characteristics.

For several years we apply now two very robust polarity-extended chromatographic separation techniques in research studies as well as in routine applications, i.e. the serial RPLC-HILIC coupling and the (polar stationary phase) SFC [2]. These two techniques provide orthogonal and complementary separations [3] and cover the full spectrum of polarity from non-polar to very polar compounds. As a consequence, the use of these techniques in environmental analysis closes the often-observed 'monitoring gap'. Both techniques offer the chance to widen the analytical scope significantly towards very polar compounds and to complement already used analytical techniques for polar and non-polar compounds such as RPLC without the requirement of two individual analysis of a sample.




Impressive examples using this chromatographic strategy presented the studies of sulfamethoxazole [4], diclofenac [5], and bisphenol A [6] degradation (in wastewater treatment) by electrochemical oxidation (in which very polar molecules like oxalic acid, fumaric acid, and/or inorganic by-products emerged).

The benefits of this 'polarity extended chromatography' have also been reported in various further matrices like red wine [7], plant metabolome [8] and house dust [9]. This gives hope for the global and fast usage of polarity extended analysis of PMTs and vPvMs monitoring in close future (especially in environmental non-target screening [10]).

References:

- [1] https://link.springer.com/journal/216/topicalCollection/AC_b22e0c8da1d84439f95dd52606862258/page/1 (latest view 06.02.2020).
- [2] Bieber S, Greco G, Grosse S, Letzel T (2017) RPLC-HILIC and SFC with Mass Spectrometry: Polarity-Extended Organic Molecule Screening in Environmental (Water) Samples. *Anal Chem* 89:7907–7914. <https://doi.org/10.1021/acs.analchem.7b00859>
- [3] Bieber S, Letzel T (2018) Orthogonal separation techniques to analyze (very) polar molecules in water samples: Assessing SFC and Reversed-Phase LC–HILIC. *LCGC Europe* 31(11): 602-606. <http://www.chromatographyonline.com/orthogonal-separation-techniques-analyze-very-polar-molecules-water-samples-assessing-sfc-and-revers>
- [4] Rajab M, Heim C, Greco G, Helmreich B, Letzel T (2013) Removal of sulfamethoxazole from wastewater treatment plant effluents by a boron-doped diamond electrode. *Int J Environmen Poll Sol* 1(3): 88-97., https://www.academia.edu/26501701/Removal_of_Sulfamethoxazole_from_Wastewater_Treatment_Plant_Effluents_by_a_Boron-doped_Diamond_Electrode
- [5] Rajab M, Greco G, Heim C, Helmreich B, Letzel T (2013) Serial coupling of RP and zwitterionic hydrophilic interaction LC-MS: Suspects screening of diclofenac transformation products by oxidation with a boron-doped diamond electrode. *J Sep Sci* 36(18): 3011-3018. <https://onlinelibrary.wiley.com/doi/10.1002/jssc.201300562>



- [6] Rajab M, Heim C, Letzel T, Drewes JE, Helmreich B (2016) Electrochemical oxidation of Bisphenol A by a boron-doped diamond electrode in different water matrices: Identification of transformation products and evaluation of inorganic by-products. *Int J Environmen Sci Tech* 2016: 13, 2539-2548. <https://link.springer.com/article/10.1007/s13762-016-1087-z>
- [7] Greco G, Grosse S, Letzel T (2013) Serial coupling of reversed-phase and zwitterionic hydrophilic interaction LC/MS for the analysis of polar and nonpolar phenols in wine. *J Sep Sci* 36:1379–1388. <https://doi.org/10.1002/jssc.201200920>
- [8] Wahman R, Grassmann J, Schroeder P, Letzel T (2019) Plant Metabolomics via RPLC-HILIC-TOF-MS. *LCGC N. America Mass Spectrom* 3: 8-15. <http://www.chromatographyonline.com/plant-metabolomic-workflows-using-reversed-phase-lc-and-hilic-esi-tof-ms>
- [9] Rostkowski P, et al. (2019) The Strength in Numbers: Comprehensive Characterization of House Dust using Complementary Mass Spectrometric Techniques, *Anal Bioanal Chem* 411 (10), 1957–1977. <https://link.springer.com/content/pdf/10.1007%2Fs00216-019-01615-6.pdf>
- [10] Schymanski E L, et al. (2015) Non-target screening with high resolution mass spectrometry: Critical review using a collaborative trial on water analysis. *Anal Bioanal Chem* 407 (21): 6237-6255. (especially see Figure S1 in Supplementary part). <https://link.springer.com/article/10.1007%2Fs00216-015-8681-7>

please cite as:

Bieber S and Letzel T (2020) White Paper – Polarity-Extended Chromatography, AFIN-TS Forum; February (1): 1-4.

AFIN-TS GmbH

Am Mittleren Moos 48

D-86167 Augsburg Germany

www.afin-ts.de

info@afin-ts.de

© Copyright 2020 – All content, especially texts, photographs and graphics are protected by copyright. All rights, including reproduction, publication, editing and translation, are reserved, AFIN-TS GmbH.

